

Monitoring Chart

QUICK TESTS™
The LEADER in User Friendly Wine Chemistry Tests

www.ACCUVIN.com

pH RECORD

Block _____

Dates	Test Results	(pH units)											
		3.0	3.1	3.2	3.3	3.4	3.5	3.6	3.7	3.8	3.9	4.0	

See Insert for more information and references. pH goals-gray.

Attach to Batch File.

COMMENTS:

pH	SO ₂ (ppm)												
	8	14	20	27	35	42	50	60	75	90	110	130	
2.9	7	8	11	13	16	21	26	32	39	49	63	79	98
3.0	8	10	13	16	21	26	32	39	49	63	79	98	123
3.1		10	13	16	21	26	32	39	49	63	79	98	123
3.2			13	16	21	26	32	39	49	63	79	98	123
3.3				16	21	26	32	39	49	63	79	98	123
3.4					20	26	32	39	49	63	79	98	123
3.5						25	32	39	49	63	79	98	123
3.6							31	39	49	63	79	98	123
3.7								39	49	63	79	98	123
3.8									49	63	79	98	123
3.9										62	79	98	123
4.0											78	98	123

Free SO₂ goals in GRAY